

Cavex Amalgam Capsules mixing times

Amalgamator		Cavex Non Gamma-2			Cavex Octight			Cavex Avalloy	
		Press	Safe	New	Press	Safe	New	Press	New
AUTOMIX	Kerr	8 - 12	5 - 7	5 - 7	8	8 - 12	8 - 12	10 - 12	10 - 12
DUOMAT (4000rpm)	Degussa	7	5 - 7	7	6	8 - 12	8	10	8
MIXOMAT	Degussa	8	5 - 7	5 - 7	6	8 - 12	8	8	8
SILAMAT (4500rpm)	Vivadent	7	5 - 7	5 - 7	6	5 - 7	5 - 7	7	7
CAPMIX	Espe	7	5 - 7	5 - 7	6	5 - 7	5 - 7	7	7
SILAMAT PLUS	Vivadent	7	5 - 7	5 - 7	6	5 - 7	5 - 7	7	7
Ultramat II	SDI	7	5 - 7	5 - 7	6	5 - 7	5 - 7	7	7
Tac 400	Linea Tac	7	5 - 7	5 - 7	6	5 - 7	5 - 7	7	7

	Safe cap	Press cap	New Capsule
Size in mm:	30,6	32,8	31,9